

Feel The WOODWAY Difference.

For The Long Run  
**WOODWAY**



Treadmills for Medicine



**WOODWAY Slat Belt Design – Superior Treadmill Technology.**

Since 1975, WOODWAY has focused on manufacturing and continuous improvement of just one product – professional treadmills.

The patented slat-belt technology, specifically designed for the human body, is clearly superior to any other treadmill technology and offers essential benefits:

**Better for your body:**

- Joint-protecting point-elastic slat-belt technology
- Smooth, slip-free acceleration from 0.1 km/h with 0.1 km/h increments
- Unparalleled running experience

**Better for the environment:**

- Significantly reduced energy consumption (low-friction ball bearing transportation system)
- Longer treadmill life
- Factory renewal program

**Better for your bottom line:**

- Considerably reduced maintenance, spare-parts and repair costs (no replacement of deck or running surface required)
- Significantly reduced operational costs
- Usage possible for patients of all sizes (maximum patient weight up to 360 kg)

Full Compliance with the Legal Requirements for Medical Products in the European Economic Area (CE marking)  
WOODWAY PPS Series treadmills comply with the Medical Devices Directive (directive 93/42/EEC, last amended by 2007/47/EC)



**Usable Running Surface**  
(W) 43/55/70 cm x (L) 157 cm,  
60 individual rubber slats

**Speed**  
0 – 20 km/h

**Elevation**  
0 – 20%

**User Weight Capacity**  
180 kg Running  
360 kg Walking

**Drive System**  
114 precision ball bearings,  
toothed belt and drive wheels

**Motor**  
2 hp continuous (5 hp peak)  
brushless servo, 1,5 kW

**Weight**  
210 – 260 kg

**Voltage Supply**  
220 – 240V 50/60Hz 16A Fuse

The WOODWAY PPS Med series and its proven history in medical application has become the first choice among hospitals, clinics and rehabilitation centers around the globe. Through continuous technology advancements and high-end material selection, including an uninterrupted medical grade handrail, WOODWAY provides each patient with a safe grip in any situation.

The WOODWAY PPS Med treadmill establishes the base for the WOODWAY PPS series product line which combines multiple customizable options across 43 cm, 55 cm or 70 cm running surface widths.





The WOODWAY PPS Plus treadmill is the ideal device for patients who require product flexibility for therapy. The fully adjustable gas assist handrails are infinitely customizable in height and width to provide a security level unmatched by more traditional stationary parallel handrails.

The generous number of features included in the WOODWAY PPS Plus series offers maximum variability providing solid comfort and accurate support for patient and therapist alike.

**Usable Running Surface**  
 (W) 43/55/70 cm x (L) 157 cm,  
 60 individual rubber slats

**Speed**  
 0 – 24 km/h

**Elevation**  
 0 – 25%

**User Weight Capacity**  
 180 kg Running  
 360 kg Walking

**Drive System**  
 114 precision ball bearings,  
 toothed belt and drive wheels

**Motor**  
 2 hp continuous (5 hp peak)  
 brushless servo, 1,5 kW

**Weight**  
 210 – 260 kg

**Voltage Supply**  
 220 – 240V 50/60Hz 16A Fuse





**Usable Running Surface**  
(W) 43/55/70 cm x (L) 157 cm,  
60 individual rubber slats

**Speed**  
0 – 20 km/h

**Elevation**  
0 – 20% (optional)

**User Weight Capacity**  
180 kg Running  
360 kg Walking

**Drive System**  
114 precision ball bearings,  
toothed belt and drive wheels

**Motor**  
2 hp continuous (5 hp peak)  
brushless servo, 1,5 kW

**Weight**  
210 – 260 kg

**Voltage Supply**  
220 – 240V 50/60Hz 16A Fuse

Due to its open handrail design, the WOODWAY PPS Ortho treadmill offers versatile applications for gait analysis. Therapist operation via serial interface and PC connection allows patients to focus on specific instruction related to treatment or diagnostics without the distraction of a monitor.

The deep black finish of the WOODWAY PPS Ortho series prevents potentially disruptive light reflections and is perfect for video applications. A reverse function provides therapists the option to record video from four different visual planes with just two cameras.





**Usable Running Surface**  
(W) 55/70 cm x (L) 157 cm,  
60 individual rubber slats

**Speed**  
0 – 20 km/h

**Elevation**  
0 – 20%

**User Weight Capacity LokoStation**  
160 kg

**User Weight Treadmill**  
180 kg Running  
360 kg Walking

**Drive System**  
114 precision ball bearings,  
toothed belt and drive wheels

**Motor**  
2 hp continuous (5 hp peak)  
brushless servo, 1,5 kW

**Weight**  
680 – 710 kg (including treadmill)

**Voltage Supply**  
220 – 240V 50/60Hz 16A Fuse

Specifically designed for training and rehabilitating gait impairments, the WOODWAY LokoStation provides effective applications for hemiplegic, paraplegic and recovering stroke patients in a safe and controlled atmosphere.

The WOODWAY LokoStation's practical design offers therapists economically feasible access to quick and convenient body weight support (BWS) including ergonomically adjustable seating for manual treatment assistance.

LokoStation Features:

- 1- or 2-point variable suspension
- Symmetric or asymmetric unloading
- Static or dynamic guidance
- Infinitely variable body weight support
- Comfortable weight suspension



## The WOODWAY Difference

In 1974 Willi Schöenberger envisioned a treadmill that was specifically designed to meet human needs - one in which the feel for the runner would emulate a wooded forest path. After intensive research with the University of Sports Science in Cologne, Germany, WOODWAY developed and patented a revolutionary machine which significantly protected the major joints used when running or walking while not interfering with natural biomechanics.

WOODWAY GmbH was founded in Weil am Rhein, Germany in 1975 with the goal to become the world's finest treadmill manufacturer. Currently, with design and manufacturing facilities in the United States, Germany and Japan, WOODWAY is proud to continue this tradition as the world's largest specialized manufacturer of sports, medical and human performance treadmills incorporating one-of-a kind slat-belt technology. Leading pro teams and coaches around the globe rely on WOODWAY treadmills to offer athletes a safe, accurate and effective workout on a machine designed for the serious rigors of professional use.

### MISSION

WOODWAY is a global organization built on the efforts and commitment of dedicated professionals. Our purpose is to manufacture products and provide service that exceeds our customers' expectations through continuous improvements in quality and design. Our focus is to building lasting customer relationships that will endure For The Long Run®.

### VISION

By viewing our customers as partners with a focus on their goals, we provide unique and superior solutions that will help our partners achieve lasting success For The Long Run®.



Authorized WOODWAY Distributor:

10/2013-V1.0 e