



## HEALING ONE STEP AT A TIME

### REHABILITATION TREADMILL

From rehabilitation to orthopedics and physical therapy, the PPS Med treadmill provides comfort and assistance to patients who require walking support. With an absolute zero start speed in increments of 0.1 km/h a patient can feel assured that they are in control and can set their own pace that best suits their comfort levels.



*\*Various widths available*

For The Long Run®  
**WOODWAY®**

WOODWAY EUROPEAN HEADQUARTERS  
Steinackerstr. 20, D-79576 Weil am Rhein, Germany  
PHONE +49 - 7621 - 940 999 - 0 | FAX +49 - 7621 - 940 999 - 40  
**WOODWAY.DE**

WOODWAY WORLD HEADQUARTERS  
W229 N591 Foster Court Waukesha, WI 53186, U.S.A.  
PHONE (800) 966 - 3929 | FAX (262) 522 - 6235  
**WOODWAY.COM**

# PPS MED



- Standard Feature
- Optional Feature

|                           |   | PPS MED  |
|---------------------------|---|--|
| USER INTERFACE            | Coded/Non-Coded Polar Heart Rate Pick-up  | ●  |
|                           | ANT+ Heart Rate Pick-up   | N/A  |
|                           | USB Power Supply  | N/A  |
|                           | Touch Grip Heart Rate Pick-up   | N/A  |
|                           | Data Monitor Display  | ●  |
|                           | Woodway User System (WUS)   | ○  |
|                           | RS 232 Serial Interface   | ●  |
|                           | Wired Remote  | ●  |
|                           | SPECIFICATIONS  | Running Surface Dimensions   |
| User Weight Capacity      |   | Run: 180 kg / Walk: 360 kg   |
| Belt Type                 |   | 60 Individual Slats  |
| Drive System              |   | 116 Precision Ball Bearings with 10 Roller Guides (4 mm lateral tolerance) |
| Running Surface           |   | Vulcanized Rubber (38-43 shore hardness)                                   |
| Drive Motor               |   | 2 hp Continuous (5 hp peak) Brushless Servo                                |
| Unit Weight               |   | 210/230/260 kg   |
| Width                     |   | 94/106/121 cm  |
| Length                    |   | 184 cm   |
| Height                    |   | 150 cm   |
| 0-20% Incline             |   | ●  |
| 0-25% Incline             |   | ○  |
| 0-20 km/h                 |   | ●  |
| 0-24 km/h                 |   | ○  |
| Reverse                   |   | ○  |
| 115 Vac 20 Amp            |   | N/A  |
| 220-240V 50/60 Hz 16A (C) |   | ●  |
| Warranty                  | 7 Year Frame<br>5 Year Function of the Running Belt<br>3 Year Drive, Motor and non-friction components<br>1 Year Full Warranty<br><br>Important:<br>The above warranty is only valid in conjunction with the performance of a WOODWAY-certified annual maintenance program. |  |

Specifications subject to change without written notice. Warranties may vary - Please contact your dealer for details.



## WOODWAY SLAT BELT RUNNING SURFACE

- Individual slats with vulcanized rubber surface
- Zero stretch & Zero slip toothed, wire-reinforced lateral belts
- Near frictionless precision ball bearing rail guide

## SAFETY & RELIABILITY

- Emergency stop magnet and safety lanyard
- Hospital grade circuitry
- Polar monitor circuitry – integrated heart rate monitoring
- Advanced Slat Belt shock absorption
- Life expectancy of running surface 150,000 to 200,000 miles
- Extended warranty and factory renewal program



## CONVENIENT USER CONSOLE

- LED screen with readouts monitoring speed, incline, distance, time, and heart rate
- Hand-held remote control (speed, elevation & stop)

Full Compliance with the Legal Requirements for Medical Products in the European Economic Area (CE marking)

- WOODWAY PPS Series treadmills comply with the following:
- Medical Devices Directive (directive 93/42/EEC)\*
  - Medical Standard EN 60601-1, 3rd edition
  - Medical Standard UL 60601-1 (USA)
  - Medical Standard UL CSA C22.2 No 601.1 (Canada)

**FOR ADDITIONAL INFORMATION ABOUT THIS PRODUCT  
VISIT [WWW.WOODWAY.COM](http://WWW.WOODWAY.COM) OR CALL US TO TALK TO A REPRESENTATIVE.**